

FUMUS GO

Self-sealing filtration kit



DATASHEET

ver. 2.0 2025

MEDICAL
SCANDINAVIA

Self-Sealing filtration

DESCRIPTION

Self-sealing filtration kit, right angled, medium volume for protection is used in combination with Medical Scandinavia A/S surgical smoke evacuation unit FUMUS GO. The Self-sealing filtration kit is used to avoid unintended events such as suction of patient tissue and blood into the FUMUS GO system.

IMPORTANT

Self-sealing filtration kit item no. 31001 should be changed between each operation, to prevent cross contamination between patients. Medical Scandinavia only guarantees FUMUS GO when using the self-sealing filtration kit.

PACKAGE	
Item no.: 31001 Self-sealing filtration kit Package: 1 pcs / 100 pcs. box / 1000 pcs. box	

Note: The self-sealing filtration is not sterile