

Product Description – Smoke exhaust filter

Compatibility	ISO Connections: Patient: 22M/15F Machine side: 15M/22F
Filtration Media	Synthetic combination pleat element
Viral Filtration Efficiency	>99.999%
Bacterial Filtration Efficiency	>99.999%
Moisture Loss	<18.8 mg/L after 24 hr with 500 ml Tidal Volume
Pressure drop hPa / mbar	@ 30L/min flow: 1.0 @ 60L/min flow: 2.3 @ 90L/min flow: 3.9

