IncuArch

ATOM Transport Incubator



Enables quick transportation and high quality care during emergencies.



Easy-to-observe hood



The size of the lighting lamp has been reduced in order to eliminate obstacles to sight. Unevenness of the hood surface has been significantly reduced to make it easier to give a blanket bath.

Constant monitoring of vital signs




The operation panel is placed at the side of the infant's head for easy observation of vital signs even during treatment.


Large operative field, accessible from both sides



The body size has been reduced while maintaining the same operative field range as that of conventional models. Snap-open access ports are placed at the head and foot sides of the hood to enable easy access from both ends.




Compact design



The weight of the main body is approx. 25kg.
Compact design.

*Weight when one battery is installed.




Securely holds the infant's body



An infant can be secured in a particular position by Dispo safety band (non-woven fabric) and Safety support mattress.

*Safety support mattress is optional.

Considers vibration



A structure that takes vibration into account has been adopted.



Products

Main units

[Code No.] 62101 IncuArch

Add-on units

[Code No.] 62110 SpO_2 unit (Masimo)
[Code No.] 62111 SpO_2 unit (Nellcor)
[Code No.] 62112 O_2 monitor
[Code No.] 22048 Lithium-ion battery pack

Specifications

Power source	Operating voltage range: AC100-240V±10% Power consumption: 220VA Frequency: 50/60Hz
Classification	Type of protection: Class I, Internally powered Degree of protection: Type BF - Applied part
Temp setting range	23.0 - 38.0°C (in 0.1°C increments)
Alarm	Excessive temperature, Set temperature, Fan, Skin temperature probe, System
Dimensions/Weight	Main body: W975 x D475 x H435mm Mattress: W620 x D330 x H30mm Weight: approx. 25kg (main unit with one battery)
Accessories	Dust cover.....1 Dispo safety band.....3 Guard.....2 Pneumoclean (Electrostatic air filter).....1 Skin Temp Probe, 5mm O.D. (Y) Cable 90cm For V-707.....1 Operation Manual.....1 Access port covers (for iris ports).....4 Quick Manual.....1 Power cord, 3m.....1



Please note that we may not be able to sell our products in regions where the transportation of lithium-ion batteries is difficult.
If you have any questions, please contact us.

Options & Supplies



[Code No.] 22043
HL Stand for IncuArch



[Code No.] 22044
Cabinet Stand for IncuArch

[Code No.] 22045 Bumper for Cabinet Stand
[Code No.] 22046 Safety support mattress
[Code No.] 22047 Dispo safety band
[Code No.] 22049 Multi-purpose Pole for IncuArch
[Code No.] 22050 Multi-purpose Pole for HL Stand
[Code No.] 22051 Crank Type Multi-purpose Pole for Cabinet Stand
[Code No.] 22052 Multi-purpose platform
[Code No.] 22053 Multi-purpose pole for multi-purpose platform
[Code No.] 22054 Monitor table for multi-purpose platform
[Code No.] 22055 Cylinder holder for 150L cylinder (Pole mounting type)

ATOM MEDICAL CORPORATION

https://atomed-global.com

Atom Medical Corporation maintains a continuing product improvement program, and therefore the equipment actually purchased and received may differ in specifications and/or configuration from that shown in the catalogue. And please note that some products may not be available in all regions.

Eng Eng-ne80-04 9.2024.Printed